

Hao blong mekem wan rapid antigen test (RAT) long nus

KEEP IT
COVID
safe


I gat ol difren kaen RAT kit we i aveilobol. Mekem sua se yu folem ol instraksen we i stap insaed long test kit.

Hemia samfala jenerol gaedlaen blong folem:


1. Ridim gud ol instraksen we i stap insaed long test kit bifo yu stat. Samfala test instraksen i gat QR kod we yu save skanem blong wajem wan video abaot test ia.



2. Wasem han blong yu.




3. Openem paket blong hem.



4. Pusum swob ia i go dip insaed long nus blong yu samples 2 sentimita (cm) mo waepem i go raon i fas long nus insaed. Mekem hemia 5 taem o blong 15 seken insaed long wanwan saed nus.



5. Putum swob ia long tiub blong wota we oli provaedem. Miksimap gud insaed. Yu save skwisim swob o brekem stik blong hem - ridim ol instraksen blong test.




6. Sarem tiub wetem kap blong hem.



7. Skwisim smol wota blong tiub i go long sampol windo blong test (hemia smol hol). Instraksen long test kit bae i talem hamas drop wota blong skwisim i go. No muvum o leftemap test kit.



8. Onem taema blong klok. Wet kasem stret taem we i talem long ol instraksen blong test kit ia blong yu mas wet. No muvum o leftemap sampol windo blong test ia.



9. Ridim risal blong hem (lukluk long raet saed pej ia).

Hao blong ridim risal blong yu



'C' Soem (Kontrol) (wan laen long ples ia hemi some se test hemi bin wok)

'T' Test (sapos i nogat laen long ples ia = negatif, sapos i gat laen long ples ia = positif)




NEGATIF

I minim se maet yu nogat COVID. Sapos yu gat eni simtom blong hem, long nekis dei bae yu mekem wan RAT bakegen o go mekem wan PCR test.



I NO WOK

Mekem wan moa test (mo ripitim ol step 2 kasem 9). Yu gat eni problem? Go mekem wan PCR test.



POSITIF

Yu gat COVID. Sapos test risal hemi soem wan laen long 'T', nomata se laen ia i smol o wik be i minim se hemi wan positif risal.

Sapos risal blong test blong yu hemi positif, yu mas aesoletem yu wan kwiktaem mo talemaot long Pablik Helt.

Skanem kod ia blong letem olgeta i save long positif test risal blong yu long intanet, o kolem olgeta long Pablik Helt Hotlaen 1800 671 738.

